

2015-08-17

how to view images and quantify signal from IVIS 20150817

View images and Quantify Signal from IVIS

Computers **INVIVO-ANALYSIS** and **analiza-4** are normally never shut down and are accessible on the Network from all computers in the building under **Microsoft Windows Network → Medicine**. To access, in Windows go to Start → Run and enter **\\INVIVO-ANALYSIS \\IVIS_Data**. Please back up your data from **INVIVO-ANALYSIS** or **analiza-4** to your computer or to your Dropbox (or similar) account.

Data from **IVIS** are stored on **INVIVO-ANALYSIS** in the following location:

C:\ IVIS_Data \[IVIS year]\mm-yyyy\[PI]

Individual experiment folders saved on **IVIS** contain a PNG file (which can be opened with common imaging or system software programs) along with TXT files of experiment information and settings.

Data set folder and PNG file are named according to PI initials, acquisition date and time:
PIyyyyymmddhhmmss

Always use locally stored file copies when you perform analysis.

Use **Living Image 3.2**, an offline software version.

When prompted to **Select User** use your PI's initials.



Visualization of experiments

File → Browse

You may choose a folder which contains a series of data set (experiment) folders or one data set folder.

In the ensuing table you may add to the list experiments performed on a different session by checking **Add to List** and pressing **Browse...**



From the list, using <ctrl> or <shift>, choose individual experiments.


Press **Load** to open an individual window for each experiment or press **Load as Group** to open a sequence of experiments in one window.

Load as Group Visualize a sequence of experiments

Press **Load as group** to open a group of experiments in one window.


(**Bug:** Part of the image window may be solid black. Resizing slightly by dragging the lower right corner corrects this).

Change Units: to **Photons**



Unify the pseudocolor scale: In the **Tool Palette** that appears, in **Image Adjust**, in **Color Scale Limits** uncheck **Individual**.

If there are great differences between the images consider checking **Logarithmic Scale**.



File → **Save as...**

The sequence is saved as a separate folder with all the sub data sets copied inside.


The image can be exported as TIFF or JPEG by pressing

The images can be viewed individually and used for quantitation by pressing

Load Visualize individual experiments

Use **Load**, or, alternatively, use **File** → **Open** and inside each experiment folder choose any file to open individual experiments side by side for visualization and/or quantitation.

In the image window change **Units:** to **Photons**.




Unify the pseudocolor scale:

Press **Info** to see the experiment parameters.


You can edit info by **Edit** → **Image Labels**

In **Tool Palette**, **Image Adjust**, **Color Scale**, select **Manual**.



From the **Color Bar** legends of the images (legend below the pseudocolor bar to the right of each image) note the lowest minimum and the highest maximum of the group.

In **Tool Palette**, **Image Adjust**, **Color Scale**, copy-paste these values at the **Min:** and **Max:** sliders of each image.



Quantitation of Regions of Interest (ROI):

Luminescent signal:


When the active window is a single experiment, the **tool palette** is extended.



Manually defined ROIs:

In **Tool Palette**, **ROI Tools**, choose circle, rectangle or contour to define one ROI at the first image (if you choose a number greater than 1 you will not be able to resize equally).

Resize or change proportions by dragging ROI handles (↖↘), reposition by panning (↕), rotate by <right-click> on the ROI, choosing Rotate, and dragging ROI handles (↻).



<Right-click> Duplicate ROI.

Reposition/manipulate duplicates.


Copy ROIs including their position by <right-click> Copy all ROIs.

In the other experiment windows <right-click> Paste ROIs.

Automatically defined ROIs:

To automatically draw ROIs detected by the software, in **Tool Palette**, **ROI Tools**, in one of the circle/square/contour choose Auto All.

To automatically draw one ROI at a user-specified location choose Auto 1. Use the ring that appears to position it.



Press **Create**.

For automatically generated ROIs define the minimum % of peak pixel intensity to be included automatically by the **Threshold%** slider in **ROI Tools**.

Rule of thumb: Use 20% threshold.

The sum of **animal non luminescent background** and **system background** signals should be subtracted from the **signal**.

To this end, the ROI of interest should be duplicated twice:


- 1) on a non luminescent place on the animal and
- 2) on a black surface in the image.

Positions of automatically created ROIs are locked.

<Right-click> Unlock positions to unlock and move.


View → ROI measurements (or click on the  icon in **ROI Tools** to obtain a table of measurements).

(You may choose to order the list by clicking on the ROI column to facilitate comparison between different experiments on the same animal.)



Average Radiance, in photons/second/cm²/steradian is a commonly used measurement of interest in luminescence.

The table can be exported by pressing Export and choosing CSV which is opened by EXCEL.



Fluorescence:

Browse to open images to create a sequence (see “visualization of experiments” above for details).


File → Save as...

The sequence is saved as a separate folder which includes all the sub data sets.




A tool tip with the excitation and emission filters used for acquisition appears when you point to the different images.

With the saved sequence file open we can use **Tools → Image Math** on Plyyyymmddhhmmss_SEQ.



We can e.g., subtract the signal acquired with the **DsRed ex. DsRed em. Bkgd.** from the one acquired with **DsRed ex. DsRed em.**



The subtracted image can be analyzed and quantified like luminescence experiments.

Online manual of Living Image 4.0

http://www.perkinelmer.com/CMSResources/Images/44-135288BRO_LIV-BR-01-5832.pdf