#### Bioimaging Center Biomedical Core Facility Ruth & Bruce Rappaport Faculty of Medicine Technion - Israel Institute of Technology



מרכז דימות ומיקרוסקופיה מרכז תשתיות ביורפואי הפקולטה לרפואה ע"ש רות וברוך רפפורט

הפקולטה לרפואה ע ש רות וברוך רפפורנ הטכניון - מכון טכנולוגי לישראל

Instructions TL1 2018-02-21.docx 2018-02-21

# **Time Lapse 1 Instructions**

Table	of	contents
-------	----	----------

Introduction	2
System components	2
Hardware:	
Objectives:	
Filter turret:	
Touchpad	
Illumination:	
Start up	
Sample mounting and viewing	5
Software	7
Acquisition	7
Auto save Tool	
Acquisition Mode Tool	9
Imaging Setup Tool	
Channels Tool	
Snapshot	
Dimensions tab	
Graphics tab	
Display tab	
Multidimensional acquisition	
Time lapse experiments	
Time Series tool	
Experiment with positions	
Tiles tool	
Sample Carrier	
Advanced setup	
Single Positions	
Position arrays	
Parameters for Time Series with Positions experiment	
Focus strategy tool	
Software Autofocus tool	
Tiling	
Local focus surface	20
Z-stack	22
Processing	23
Time concatenation	23
Tile stitching	23
Shading	23
Image Export	
Export cropped ROI	
Split scenes (positions)	25
Movie export	25
Batch	
Batch method tool	25
Shutdown	

## Introduction

This is an inverted widefield system wholly enclosed in a polycarbonate incubator system with controlled temperature , humidity and an active gas mixer used mainly for short or long term time lapse experiments.



## System components

Microscope: Inverted Zeiss Observer Z1 with motorized stage inside a BOX incubator, specifically designed to encompass the specific microscope BRICK gas mixer that delivers and monitors N<sub>2</sub>, Air, CO<sub>2</sub> CUBE temperature controller Humidifier Mini chamber that delivers humidified mixed gas to the microplate/dish N<sub>2</sub>, Air, CO<sub>2</sub> cylinders (outside of the room) Computer, monitor, UPS Camera: AxioCam HR R3



Hardware:						
Objectives:						
Objective	Magnification	NA	contrast	Dish bottom	Working	Immersion
				or slide coverslip	distance (mm)	
EC Plan-Neofluar	x5	0.16	Phase 1	plastic/glass	18.5	air
EC Plan-Neofluar	x10	0.3	Phase 1	plastic/glass	0.3	air
Plan-Neofluar	x20	0.4	Phase 2	plastic/glass	0.4	air
Apochromat	x20	0.8	Phase2 & DICII	Glass	0.55	air
Apochromat	x40	0.95	Phase 3	Glass	0.25	air
Plan-Apochromat (optional)	x63	1.4	Phase 3	Glass	0.19	oil

## Filter turret:

Reflector Turret	Excitation	Beamsplitter	Emission	Suitable fluorophores
Analyzer module DIC				
Zeiss Filter set 50	BP 640/30	FT 660	BP 690/50	Cy5, DRAQ5
Zeiss Filter 20HE	BP 546/12	FT 560	BP 607/80	PI, Cy3, Alexa555, Rhodamine
Zeiss Filter set 38HE	BP 470/40	FT 495	BP 525/50	Cy2, GFP, Alexa 488
Zeiss Filter set 49	G 365	FT 395	BP 445/50	DAPI, Hoechst
Optovar x1.6 (optional)				

## Touchpad

The microscope can be controlled by the touchpad: objectives, reflector turret filters, transmitted or reflected illumination, light path. A coaxial focus wheel is attached.

## Illumination:

Transmitted:HalogenEpifluorescence:X-cite metalhalide



## Start up

- 1. Check both gas tanks' pressure; don't start a new live cell experiment if pressure is below 30 bar. If in doubt, please ask for assistance from BCF staff.
- 2. Turn each gas tank's black spigot CCW to open-horizontal (please don't move any valve).



3. Turn on **BRICK**, check parameters (1).



- 4. Make sure **CUBE** temperature is set to 37C°.
- 5. Check the water level in the humidifier; to refill please ask for assistance from BCF staff.



Turn on X-cite lamp if you need fluorescence (2).
 Check intensity setting, set it to a low setting for imaging live cells.
 Do not turn off within half an hour of turning it on.



7. Turn on microscope controller (**3**).





8. Turn on microscope (ON\OFF button on left side of microscope behind a little door) (4).



- 9. Turn on computer and monitor
  - Username: multilabs
  - Password: 123456
- 10. Log in to your BookItLab account and activate your reservation.
  - The ZEN 2.3 software will subsequently open.
- 11. In the software window allow stage calibration.

#### Calibrate Now

Please note, the stage will move automatically; make sure there are no objects interfering in its path. If calibration dialog does not appear, close the software and restart the computer.

## Sample mounting and viewing

Insert your plate and make sure both stage insert and plate sit firmly and do not wobble.

Cover with the mini chamber lid.



Manually choose the condenser position according to the objective (see Objectives Table).





## Choose light path $\rightarrow$ eyepiece/VIS

Home	tives flector path	
	Sideport L	Load position
Control		TL Illumination
Automatic	100%	Off On
		RL Illumination
XYZ	0%	Off On
Display		

Focus on your specimen.

If you acquire transmitted light you will have to set up Kohler illumination:

- **1.** Bring the sample into focus.
- 2. Close the field diaphragm until you can see at least one edge.
- 3. Adjust the condenser height until the edges of the diaphragm image are crisp.
- 4. Center the diaphragm image using the two centering screws.
- 5. Open the field diaphragm, just until the image fills the field of view.















## Software

Zen Blue is organized such that on the left, below the menus there are four **main tool tabs** with their respective tools below



In the middle there is the **image container** with tabs for the different open images, **Live** view and Tile/Position -Advanced Setup View with more tabs below

On the right there is the open file list and a Devices section to change/move Microscope objectives, Stage and Focus.





Load a saved configuration which includes the desired channels and acquisition settings, such as multidimensional and saving options.

- \*-

Revital BF Cy2



## Alternatively, use Smart setup to define channels:

Click on 🗡 Smart Setup

To add a channel click on the [+] button at the **Configure your experiment** section.



## The Add dye or contrasting method list appears

Double-click on the desired fluorophore or TL Brightfield under Contrast methods for transmitted.

Repeat for further channels.

Click ok

It is strongly recommended to use **Auto save** for time lapse experiments.

In the experiment options below the start buttons check auto save

Experiment Regions	Experiment Designer
Auto Save	Automated Export
Automation	

The Auto Save tool will appear in the Applications section below



Choose the folder your files will be saved in and give a name.

This will be used as a prefix in case more than one files are created before altering the **Name** field and automatic numbering will be appended.

Save your files in:

D:\Users data\PI name\User name

Checking Automatic sub-folder will create a folder YYYY-MM-DD where your files will be saved.



Check settings:

**Acquisition Mode Tool** 

## In Camera section:

Usually no binning is required (1x1) and maximum ROI is acquired (Resolution/Size 1388x1040)



#### Gain should be set to 1; no post processing actions

<ul> <li>Gain</li> <li>Analog Gain</li> <li>1x</li> </ul>		
<ul> <li>Post Processing</li> <li>Black Reference</li> </ul>	Define	
Enable Noise Filt Threshold	er	1.5
Enable Unsharp I Strength Radius Color Mode	Mask 0 Luminance	1.0
☐ Auto Contrast Contrast Tolerance ✓ Clip to Valid Bits	-0	1.0

## For accuracy, Slow live speed and High accuracy readout speed Only $\checkmark$ NIR mode if you are acquiring infra-red channel (Cy5)

Mode		
Live Speed Slow	•	
NIR Mode		
• Trigger Control		
Model Specific		Reset
General		
	Original	•
	High accuracy (12.5 MHz)	•



## **Imaging Setup Tool**

After setting each channel as specified above there is no need to touch the Imaging Setup configuration.



#### **Channels Tool**

Press Live in the start buttons



This will open a Live image tab in the center section.

Use the tabs below the image:

Dimensions to control zoom and choose to display in **range indicator** and

Display to change brightness/contrast.

For each channel you will have to:

•  $\sqrt{}$  the channel and click to select it (like "AF647" in the image below)

Bright	Ref.	•
✓ AF647		•
Set the exposure time using the slider		
Time 204.0	<b>•</b>	ns 🔻

For live experiments: Shift 70% for transmitted light, 30% for fluorescence channel
 Shift 70%



## Snapshot

To acquire a single snapshot check the pertinent channels in the Channels tool

🔽 Bright	Ref. 🗖 🏲
AF647	<b>—</b> •
AF555	🗖 🔻
AF488	<b>•</b> •
AF405	<b>□</b> ▼
And press Snap in	the start buttons



If you do not **auto save**, save the file using the **Images and documents** tool at the right panel.

Save file as native **.czi** which contains all the experiment parameters and could be requested by any reviewer at time of publication. Always backup and keep the native file.

A file is not saved if it has an asterisk after its name.



To change a saved file's name right click on the file's name at the top of the tab and choose Rename

Ne			
INC	Save	Ctrl+S	
	Save As CZI	Ctrl+Shift+S	
(	Rename		
	Delete		
	Close	Ctrl+F4	
	Close All	Ctrl+Alt+W	
	Close All But This		
	Open Containing Folder		

#### Tabs below the image:

#### **Dimensions tab**

Use the dimensions tab to zoom (alternatively you can zoom with the mouse wheel), move between different time-points and z-levels in case of multidimensional image.

 $\checkmark$  Range indicator to monitor saturation (avoid red pixels which are saturated)

Dimensions	Player	Graphics	Custon	n Graphic	s				
Time 1	0—	,				, 4	1		
Zoom [	¢ 10	0% Q	Ð	-0-		36	% 🗘	🖌 Auto F	
Tools	<b>₿</b>	گ⊡ 🖑	ø		lavigator		✓ In	terpolatio	'n
Channels	Ch2	Ch4	ChD						
	Single	Channel	Ű	ange Ind	licator	Q	uick Co	lor Setup	



#### **Graphics tab**

Use the graphics tab to add annotations, e.g., scale bar or relative time.

Press format or right click on the annotation to change format



#### **Display tab**

Change image brightness and contrast by selecting a channel and moving the line defining the limits of gray levels. Note that in order to compare images it is absolutely imperative to change the parameters to exactly the same levels, i.e., the values of Black, Gamma, White should be the same.



### **Multidimensional acquisition**

#### **Time lapse experiments**

 $\checkmark$  Time series in the multidimensional option section



The Time Series tool will appear



#### **Experiment with positions**

 $\checkmark$  **Tiles** in the multidimensional option section.



The Tiles tool will appear.



### Tiles tool

### Sample Carrier

First select your sample carrier and calibrate it



Go through the wizard:

1/7 Select Low Magnification Objective

Select x5 or x10 objective and press Next

#### Set transmitted illumination and press Next

Stage will have been calibrated at startup. Press Next (below)



For multiwell choose Search well edges (7 points)

#### Press Next



Use Move stage to go to the approximate position of the A1 well's leftmost point.

You may need to use the eyepieces to locate the relevant point and then return to **Live** to finetune. Click **Set position**, repeat for next three points, click **Next** 

6/7 Define Upper Right Well (A6)
Repeat for upper right well
Repeat for lower right well
6/6 Define Lower Right Well (D6)
Press Finish
Continue with position setup below
Return to Tile setup here.

Delete any old tiles or positions by selecting and pressing



## **Advanced setup**

Use **Advanced setup** in **Tiles** tool to set up positions for the experiment



This opens a new tab in the image container (in the middle of the screen) with the **live view** in context within the vessel and with the ability to zoom in and out with the mouse wheel.

You can set and remove positions using an array of tabs below the Advanced setup View



To move between wells, in Carrier tab double click on a well. Stage moves to the center of selected well.



## **Single Positions**

To manually define positions in **Tiles** tool choose **Single positions**.



Use Live to locate each position and press

Alternatively use the + tool in the position setup tab, "setup by location"

(the arrow tool will select and move an existing position)



### **Position arrays**

For randomly set positions in Tiles tool choose Position arrays.



In the Advanced Setup View tabs in Carrier tab (center bottom) choose wells to apply the set up.



In position setup tab choose method:

(1) Position setup  $\rightarrow$  Array  $\rightarrow$  Contour Square/Ellipse  $\rightarrow$  Number e.g., "4"

In Advanced Setup View, in one well draw a square/ellipse and 4 positions will be distributed randomly within the specified shape.





(2) Position setup  $\rightarrow$  Array  $\rightarrow$  Contour Square/Ellipse  $\rightarrow$  Grid  $\rightarrow$  Overlap e.g., "-50%"

In Advanced Setup View, in one well draw a square/ellipse and as many positions as can fit with the indicated overlap will be created.

(3) Position setup  $\rightarrow$  Carrier  $\rightarrow$  Grid  $\rightarrow$  e.g., Columns "2", Rows "2", Overlap "-20%"  $\rightarrow$  +Create



(4) Position setup  $\rightarrow$  Carrier  $\rightarrow$  Number  $\rightarrow$  e.g., "4"  $\rightarrow$  +Create These will be randomly set without a pattern

In the Tiles tool Positions section you will see the array such that for each well A1, A2... there are 4 positions P1, P2, P3, P4 set in a regular grid pattern.

Single Positions			Position Arrays				
		â	Posit	tions	Size (µm)		
	A1	OÍ	4	1	15360.4 x 1	5360.4	
$\checkmark$	A2	0 🔒		1	15360.4 x 1	5360.4	
$\mathbf{\nabla}$	A3	0 🔒		1	15360.4 x 1	5360.4	
$\mathbf{\nabla}$	A4	0 🔒		1	L5360.4 x 1	5360.4	
~	^ m	1				*	-
		ected Ari					
	Name	Categ	ory	X (µm)	Y (µm)	Z (µm)	
$\mathbf{\nabla}$	P1 +	Defaul		10048.6	13008.9	7384.0	
$\checkmark$	P2 +	Defaul		11123.0	13008.9	7384.0	
V	P3 +	Defaul		10048.6	13813.9	7384.0	
V	P4 +	Defaul		11123.0	13813.9	7384.0	
						*	•
	Verify Positions						

To move a SINGLE position, it will have to be selected in the list of the single positions or of the **Positions of selected array** and then in the **Position setup** tab choose **Contour ()** and manually move the selected position/array in the **Advanced Setup** View (changes from **+** to **()**. Position arrays will move en bloc.



Regardless of way of defining positions, they will need to be verified.

Press Verify Positions below the positions list.

The verify positions dialog opens.

Use AF to Verify the Remaining



Double-click first position, correct focus and press Set Z & Move to Next until all points have been verified. Alternatively you can use autofocus:

Use \_\_\_\_\_\_ as helper method, choose the first position in the list, press \_\_\_\_\_\_ and then

#### Parameters for Time Series with Positions experiment

In Tile tool  $\rightarrow$  Options section set parameters as below:

Tile overlap: 10-15%  $\rightarrow$  Travel: comb  $\rightarrow$  Positions: Sort by Y, then X  $\rightarrow$  container: meander  $\rightarrow$  Use stage speed from stage control  $\rightarrow$  backlash correction  $\rightarrow$  DO NOT stitch during acquisition  $\rightarrow$  Image pyramid during acquisition

Options
Tile Overlap 10 %
Stage Travel Optimization
Travel in Tile Regions Z Comb
✓ Tile Regions/Positions Sort by Y, then X ▼
✓ Carrier Wells/Container 🔁 Meander 💌
✓ Use Stage Speed from Stage Control
✓ Stage and Focus Backlash Correction
Split Scenes into Separate Files
Stitching During Acquisition
✓ Image Pyramid During Acquisition

### Focus strategy tool

In Focus Strategy tool right-click on a channel to set as reference (usually Brightfield)





#### Set parameters as indicated below:

Local (per Region/Position)  $\rightarrow$  Adapt Focus Surface/Z values  $\rightarrow$  with Software AF  $\rightarrow$  as additional action  $\rightarrow$  Initial Definition: By Tiles setup  $\rightarrow$  Stabilization: Expert  $\rightarrow$  Synchronized: Positions, Every 1 Position

C Focus Strategy ✓ Show All							
Use Focus Surface/Z Values Defined by Tiles Setup							
🗢 Re	ference Channel						
	Name		Offset (µm)				
/							
	AF647						
	AF555						
$\checkmark$	AF488						
	AF405						
Set	as Reference Channel						
○ For LC ✓ A with	<ul> <li>Focus Surface</li> <li>Local (per Region/Position) Global (Carrier based)</li> <li>Adapt Focus Surface/Z Values</li> <li>with Software AF          <ul> <li>As Additional Action</li> <li>Is Additional Action</li> </ul> </li> </ul>						
Initia By	al Definition of Support Points/I Tiles Setup	Positions/T	ile Regions	D			
🗢 Sta	bilization Event Repetitions and	d Frequenc	у				
	Standard Expert						
Synchronized with Image Acquisition							
	Entry Contraction Contrac						
V	✓ +++ Positions						
Repeat Every 1 Position							

### Software Autofocus tool

For long term experiments it is advisable to use autofocus after setting initial focus manually.

#### Also to check and correct after a few hours.

#### Use the settings below:

Mode: Auto  $\rightarrow$  Quality: Default  $\rightarrow$  Search: Smart  $\rightarrow$  Sampling: Default  $\rightarrow$  Relative Range  $\rightarrow$  Automatic Range



Press Start Experiment

If needed you can pause the experiment and verify position z again after checking how long you have before the next imaging cycle in the status bar below.



## Tiling

 $\checkmark$  **Tiles** in the multidimensional option section.



The Tiles tool will appear.

If your experiment contains positions or robotic movement between wells, you should calibrate the vessel as specified <u>here</u>.

Delete any old tiles or positions by selecting and pressing **mathematical below** the tiles and positions lists in the **Tiles tool**.

Several options to set an individual tile region:

Press Live

(1) In Tiles tool by contour square/circle, tiles and (e.g., "5x3") and click

Tile Regions						
			0			
		Tiles		Size	Stake	
	5		4118.2 μι	n 🌻		
	3		1878.2 µı	n 🗘		

(2) In Tiles tool by contour square/circle, Stake

Press **e** to set a mark at the current stage position in order to define a square or circular tile region.

Tile Regions						
		)				
	Tiles	Size	Stake			
Extend with t	further marker	s 🔫	Done			

Press Done to finish



(3) To define a polygonal tile region, in the Tile Region Setup tab below the Advanced Setup View (to the right) mark contour  $\rightarrow$  polygon and at the Advanced Setup View mark the limits of the polygon (right click to finish)



You can change the shape of the polygon by right-click  $\rightarrow$  edit points



You can also add more points to the yellow contour (to finetune) by clicking on it.

(4) To set tile regions in a number of wells, in the Carrier tab below the Advanced Setup View (to the left) choose wells.

In the Tile region setup tab (below, right) choose **Carrier** Choose one of the options below and click + Create





(5) To add tile regions repetitively use the Tile Region Setup tab below the Advanced Setup View (to the right) select **Predefined** and choose e.g., "Contour"  $\rightarrow$  "Square"  $\rightarrow$  "5x3". Then use the **Figure** button to plant the 5x3 square tile regions.



(6) To copy/paste a tile region right click on it at the Advanced Setup View and choose copy;



click on a new position and press paste.



(7) to copy all tiles from a well onto another well right-click, choose copy container for replication,



move to destination well, right-click paste container for replication



You can preview the tiles (with lower magnification objective for speed) to possibly correct placement and size. In Advanced setup View, in the Preview scan tab (below, left) press **Start Preview Scan** 



#### Local focus surface

To create focus surfaces with support points across the tile region mark focus positions which will be interpolated across the entire tile region:

Select a tile region in the Tile regions list in the Tiles tool

or click on a tile region in the Advanced Setup View.





## In the Focus surface (verify) section in the Tiles tool

Choose interpolation degree according to the size and unevenness of the sample



## In the Support Points tab below the Advanced Setup View:

To distribute randomly, define number of support points (usually 4 or 9, according to size and unevenness of sample) and press Distribute.



#### To manually add support points:

In the Advanced Setup View, double click where you want to center the stage.

In Support Points tab below the Advanced Setup View (below right)





In Tiles tool, in Focus Surface (verify) section click on Verify Tile Regions/Positions...

And follow the dialog that opens.



### **Options:**

Use the Options section in the Tiles tool to define parameters: 10%, Comb, Carrier meander,  $\checkmark$  use stage speed from stage control if you want to slow down (default is 100%).

Do not stitch during acquisition.



## Z-stack

Check Z-stack in the multidimensional section



## The Z-stack tool will appear.



### Click First/Last.

In Live mode adjust focus until you have reached the upper/lower plane of the Z-stack.

Click on Set First.

Adjust focus until you have reached the lower/upper plane of the Z-stack

Click on Set Last.

After choosing all channels, click on **optimal** interval to adapt to the Nyquist criterion according to channels and microscope configuration.





## Processing

## **Time concatenation**

Open two files

Choose Processing Main Tool  $\rightarrow$  Time concatenation (not in Zeiss Zen Lite)

In Input choose first file and second file

Press Apply

### **Tile stitching**

In **Processing Main tool**  $\rightarrow$  Geometric  $\rightarrow$  Tile stitching (not in Zeiss Zen Lite)

- In Input choose file
- In Parameters choose New Output

 $\checkmark$  Fuse tiles and  $\checkmark$  Correct shading (usually automatic will suffice)



Save new file

### Shading

Acquisition  $\rightarrow$  Live  $\rightarrow$  choose field of view without info and manually defocus

Acquire time series of 20 images

Start experiment

Save file

Processing  $\rightarrow$  time series  $\rightarrow$  gliding average  $\rightarrow$  averaging length  $\rightarrow$  20



input saved shading reference image

Press Apply

Processing  $\rightarrow$  Smooth  $\rightarrow$  low/pass filter  $\rightarrow$  input  $\rightarrow$  kernel size big (~15 and 15)  $\rightarrow$ 



Input Definition Set Input Automatically

After processing 
Switch to Output

Input image for stitching, input processed reference image



Bioimaging Center Biomedical Core Facility The Ruth & Bruce Rappaport Faculty of Medicine Technion-Israel Institute of Technology bcf.technion.ac.il Apply

## $\checkmark$ fuse tiles $\rightarrow\,\checkmark\,$ shading $\rightarrow\,$ reference

	Inplace		Ň	lew Output	
I Fusi	e Tiles				
M Cor	rect Shading	Referenc	e		
Press	Apply	'			

#### Image Export

This is included in the **Zen Blue** lite edition which can be <u>downloaded</u> from the Zeiss website.

ALWAYS KEEP YOUR ORIGINAL czi FILES

Processing Main Tool Tab  $\rightarrow$  Method  $\rightarrow$  Export/Import  $\rightarrow$  Image Export

In Input choose image.

In Parameters choose TIFF (lossless); JPEG is lossy and not suitable for scientific imaging.

Choose options such as exporting each channel plus the merged channel view, greyscale or pseudocolor, changes in display (<u>brightness/contrast</u> – note, should be the same values for all images), annotations.

✓ Parameters ✓ Show All
Settings
File type Tagged Image File Format (TIFF)
Compression <b>LZW</b>
Resize 1 10 100 %
<ul> <li>Original Data</li> <li>Sinit Pixel</li> <li>Apply Display Curve and Channel Color</li> <li>Burn-in Graphics</li> <li>Merged Channels Image</li> <li>Individual Channels Image</li> <li>Use channel names</li> </ul>
<ul> <li>Use Full Set of Dimensions</li> <li>Define Subset</li> </ul>
Export to C:\Users\omer\Desktop
✓ Create folder ☐ Generate xml file ☐ Generate zip file
Prefix MEF 100k cells per ml #2_Maximum intensity p

### **Export cropped ROI**

Press

In Graphics tab choose and draw a rectangle on the image.

In export parameters use define subset: subset  $\rightarrow$  Region  $\rightarrow$  Rectangle region

Use Full Set of Dimensions Define Subset					
Channels	All				
Time	All				
Region	Full				
Full					
• Full					
Ex Rectangle region					



## Split scenes (positions)

Processing $\rightarrow$	Method	: Split scenes (v	vrite files)		
Method					
Recently used					
Stitching					
Gliding Average	Gliding Average				
Time Concatenation	Time Concatenation				
Split Scenes (Write files)					
Image Export					
In Input tool of	choose m	aster file and p	ress		

#### **Movie export**

Processing  $\rightarrow$  Method: Export/Import  $\rightarrow$  Movie export





Choose which timepoints/channels/positions will be included:

e.g., All, Extract Single, Extract range: "1-4", Extract Multiple: "1,4,10"



Set output folder

Apply



#### Batch

Not all processing methods can be performed in batch.





#### To batch-export images press "Batch"



## and choose Image export from the Batch Method tool,

Press + Add to add files and select output folder



#### Select one file from the list



#### Set parameters in Parameters tool



#### Press Copy Parameters

#### Select all other files from the list



## Press Paste Parameters

#### Select all files

Use Input Folder as Output Folder Naming				Parameters Pa All Run Select	ste Parameters
S					Output N
2		D:\Users\Multilabs\Pictures\	11.18 MB	Image Export	
2/		D:\Users\Multilabs\Pictures\	11.18 MB	Image Export	
2		D:\Users\Multilabs\Pictures\	11.18 MB	Image Export	
2		D:\Users\Multilabs\Pictures\	11.19 MB	Image Export	_
		D:\Users\Multilabs\Pictures\	11.19 MB	Image Export	
		D:\Users\Multilabs\Pictures\	11.19 MB	Image Export	
~	✓ + Add	– Remove 🗑 Remove All		🕞 Load List	Save List



The same method can be applied to stitching (no option for shading reference, only automatic).



## Shutdown

- Check that all your files are saved properly.
- Close the software
- Log off from your BookItLab account.
- Copy your data from the local folder to the BCF server.



- Please do not use any form of USB flash disk to copy your files.
- Shut down the computer.
- Close the microscope controller (3).
- Close the X-cite lamp (2).
- Close the BRICK (1).
- Turn both gas tanks' black spigots CW to close (please do not move any valve).

