

DEMONSTRATED PROTOCOL

Removal of Dead Cells for Single Cell RNA Sequencing

Overview

A high percentage of non-viable cells may impact the targeted cell recovery in 10x Genomics Single Cell protocols. This Demonstrated Protocol outlines removal of dead cells from single cell suspensions. The protocol was demonstrated using peripheral blood mononuclear cells (PBMCs), dissociated tissue cells from colorectal cancer (CRC), and clear cell renal carcinoma (CCRC) patients. However, it may be used as a basis for removing dead cells from other primary cells as well as other cell lines in preparation for use in 10x Genomics Single Cell protocols.

Additional Guidance

Consult Demonstrated Protocol Cell Preparation Guide (Document CG00053) for Tips & Best Practices on handling cells and Technical Note Guidelines for Accurate Target Cell Counts using 10x Genomics Single Cell Solutions (Document CG000091) for determining accurate cell counts.

Input PBMCs for this protocol were thawed, washed, and counted as described in Demonstrated Protocol - Fresh Frozen Human Peripheral Blood Mononuclear Cells for Single Cell RNA Sequencing (Document CG00039).

Cells carry potentially hazardous pathogens. Primary tumor cells should be handled under biosafety level-2 (BSL-2) conditions. Follow material supplier recommendations and local laboratory procedures and regulations for the safe handling, storage, and disposal of biological materials.

Cell Sourcing

Dissociated tissue cells from CRC and CCRC patients were acquired from Discovery Life Sciences. The cells were thawed, washed, and counted per the Discovery Life Sciences protocol, except the steps describing the DNase I treatment of the cells were omitted. Refer to www.dls.com for specifics.

Preparation – Buffers

Buffers	Composition
Maintain at Room Temperature	
1X Binding Buffer	1X Binding Buffer Stock* in Nuclease-free Water
PBS + 0.04% BSA	


* Provided with the MACS Dead Cell Removal Kit as a 20X stock solution.

Specific Reagents & Consumables

Vendor	Item	Part Number
Thermo	Trypan Blue Stain (0.4%)	T10282
Fisher Scientific	UltraPure Bovine Serum Albumin (BSA, 50 mg/ml)	AM2616
	Countess II FL Automated Cell Counter	AMAQAF1000
	Countess II FL Automated Cell Counting Chamber Slides	C10228
	Nuclease-free Water	AM9937
Millipore-Sigma	Phosphate-Buffered Saline (PBS) with 10% Bovine Albumin (alternative to Thermo Fisher product)	SRE0036
Miltenyi	MACS Dead Cell Removal Kit	130-090-101
	MS Columns	130-042-201
	LS Columns	130-042-401
	MACS Multistand	130-042-303
	OctoMACS Separator (for use with MS Columns)	130-042-109
	QuadroMACS Separator (for use with LS Columns)	130-090-976
Bel-Art	Flowmi Cell Strainer, 40 µm	H13680-0040
Eppendorf	DNA LoBind Tubes, 2.0 ml	022431048
Corning	Phosphate-Buffered Saline (PBS) without Calcium & Magnesium	21-040-CV

Protocol Overview

Dead Cell Removal



*Resuspensions are in PBS + 0.04% BSA

Protocol

This protocol was demonstrated using sample sizes compatible with Miltenyi Biotec MS columns ($\leq 2 \times 10^8$ total cells). If using sample sizes compatible with Miltenyi Biotec LS columns ($\leq 2 \times 10^9$ total cells), consult the manufacturer's instructions.

- a. Centrifuge the cell sample at **300 rcf** for **5 min** (PBMCs) or **10 min** (dissociated tumor cells).
 - b. Remove the supernatant without disturbing pellet.
 - c. Add **100 μ l** Dead Cell Removal MicroBeads and resuspend cell pellet using a **wide-bore** pipette tip.
 - d. Incubate for **15 min** at **room temperature**.
 - e. During incubation, prepare MS column by rinsing with **500 μ l** 1X Binding Buffer. If using the LS column, consult the manufacturer's instructions for more information.
 - f. After incubation is complete, dilute the cell suspension (containing Dead Cell Removal MicroBeads) with **500 μ l** 1X Binding Buffer.
 - g. Apply cell suspension to the prepared column. The positively selected dead cells will be retained on the column while the negatively selected live cells pass through the column.
-  **DO NOT** apply the plunger supplied with the column, otherwise positively selected dead cells will be collected in the effluent.
- h. Collect the effluent (effluent 1) containing the live cell fraction in a 15-ml centrifuge tube.
 - i. Rinse the column with **2 ml** 1X Binding Buffer and collect the effluent (effluent 2). Combine effluents 1 and 2. If using the LS column, consult the manufacturer's instructions for more information.
 - j. Centrifuge cells at **300 rcf** for **5 min** (PBMCs) or **10 min** (dissociated tumor cells).
 - k. Remove the supernatant without disturbing the pellet.
 - l. Using a **wide-bore** pipette tip, add **1 ml** 1X PBS + 0.04 % BSA to each tube and gently pipette mix 5x to resuspend cell pellet.
 - m. Transfer the cell suspension to a 2-ml tube.
 - n. Centrifuge cells at **300 rcf** for **5 min** (PBMCs) or **10 min** (dissociated tumor cells).
 - o. Remove the supernatant without disturbing the pellet.
 - p. Resuspend pellet in **1 ml** 1X PBS + 0.04 %.
 - q. **Repeat** n-o for a total of two washes.
 - r. Resuspend pellet in 1X PBS + 0.04 % BSA using a **regular-bore** pipette tip to achieve a cell concentration of **700-1,200 cells/ μ l**. Gently pipette mix 10x or until cells are completely suspended.
 - s. Determine the cell concentration using a Countess II FL Automated Cell Counter. If necessary, dilute the cells with additional 1X PBS + 0.04% BSA until the target cell concentration is reached.
 - t. Once the target cell concentration is achieved, place the cells on ice.
 - u. Proceed **immediately** to the 10x Genomics Single Cell protocols.

Results

To demonstrate the efficiency of this protocol, a PBMC sample containing a high percentage of non-viable or dead cells was used. The results obtained following this protocol are outlined below. To determine percent viability, cells were stained with trypan blue and counted using a Countess II FL Automated Cell Counter.

Sample	Total Cells Before Removal	% Live Before Removal	% Live After Removal
PBMCs	1.0×10^6	49%	85%
Colorectal cancer	6.5×10^6	55%	80%
Clear cell renal carcinoma	3.7×10^6	30%	70%

© 2020 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at <http://www.10xgenomics.com/legal-notices>, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of 10x Genomics products in practicing the methods set forth herein has not been validated by 10x Genomics, and such non-validated use is NOT COVERED BY 10x GENOMICS STANDARD WARRANTY, AND 10x GENOMICS HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner 10x Genomics terms and conditions of sale for the Chromium Controller or the Chromium Single Cell Controller, consumables or software, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by 10x Genomics that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact:
support@10xgenomics.com
 10x Genomics
 6230 Stoneridge Mall Road
 Pleasanton, CA 94588 USA

