



## DNA quantity (Qubit/picogreen)

### Sample Delivery Instructions

1. Please deliver 5 $\mu$ l DNA at a concentration range of 0.5-100ng/ $\mu$ l (for Qubit) or 2-150ng/ $\mu$ l (for picogreen)
2. Samples should be delivered in 8 tubes strips or in 96 well plate.
3. An electronic sample information sheet should be filled and sent by e-mail.
4. Please give the samples a number (1, 2, etc.) written on side of the tube (same number as in the sample sheet) in case samples are delivered in 8 tubes strips. If using a 96 well plate please specify the samples position in the sample information sheet.