



RNA Quality

Sample Delivery Instructions

1. Please deliver 3 μ l RNA at a concentration range of 30-200 ng/ μ l (regular assay) or 0.5-25 ng/ μ l (high sensitivity assay).
2. Samples should be delivered in 0.2ml tubes or in 8 tubes strips.
3. Deliver samples frozen (dry ice) to assure sample integrity.
4. An electronic sample information sheet should be filled and sent by e-mail.
5. Please give the samples a number (1, 2, etc.) written on top of the tube (same number as in the sample sheet).